


► Redesigned – Experience the ease of use

pH-Fix 0-14 PT – snap cap tubes



NEW!

Improved: Simple one-hand operation

The snap cap of the tube can be **opened and closed** just **with the thumb** of the holding hand. This makes working with pH-Fix 0-14 PT comfortable and relaxing.



Improved: Relaxing easy-close principle

Due to the **long tube** design the strips do not stick out. This effectively **prevents blockades of the cap** due to jammed strips. The tube can easily be closed with a push of the thumb. This makes pH-testing convenient and relaxing.



Improved: Maximum safety concept

Geometry and material of the tube have been specially selected for maximal stability. The tube is **virtually unbreakable**. It stands firmly on its circular bottom, even when it is slightly moved. Working with these new tubes is safe.

